

Microcentrifuge



- Additional PCR strip rotor (2 x 8 x 0.2mL) and reduction adaptors for 0.2/0.4mL microtubes supplied with the product
- Electronic safety brake provides an immediate stop upon opening the lid
- 8 slot closed rotor for quiet and powerful operation
- Quick operation by simply closing and opening the lid
- Digital display of speed and timer function (CR-68X)
- Imbalance cut-off safety feature (CR-68X)
- Digital calibration function (CR-68X)

Ordering Information 2025

Cat. No.	Description	
CR-68	CAPPProndo Microcentrifuge 6.000 RPM/2.000g	152.00 €
CR-68X	CAPPProndo Microcentrifuge w/adjust. Speed & timer function 6.000 RPM/2.000g	179.00 €

Universal Centrifuge

CRC-432X 1 457.00 €
CAPPProndo Universal Centrifuge



CR-43-PL 322.00 €
2 x microplates swing rotor, RCF 2355 x g



CR-43-15F-24 243.00 €
24 x 15mL fixed-angle rotor, RCF 2852 x g



CR-43-50-4F Price on request.
4 x 50mL fixed-angle rotor, with stainless steel RCF 2852 x g



Specifications :

Speed Range	500 – 4.500 RPM/ (g) varies per rotor type
Rotors (optional)	Multiple rotor options (swing-out & fixed-angle rotors)
Remote Operation	Yes, via USB
Display	Yes, large LED Display
Timer Setting	From 1 to 999 mins or infinite mode
Programmable	Yes, up to 99 programs
Imbalance Detection	Yes
Dimension (WxDxH)	595 x 436 x 311mm
Weight	23 kg (without rotor)

More rotors available, request our catalogue.

Magnetic Stirrer

CRS-15X 158.00 €
CAPPProndo Magnetic Stirrer
1.500 RPM w / 1 Stirring bar



Magnetic Hotplate Stirrer

CRS-21H 294.00 €
CAPPProndo Magnetic Hotplate Stirrer
2.200 RPM, 320°C ceramic, 10L capacity



3D Shaker

CR3-08X 657.00 €
CAPPProndo 3D Shaker

Plate size 350x350mm
Speed range 10-120 rpm



A&E Scientific srl – bv

Rue de Lekernay, 11 / Lekkernijstraat, 11

B-7850 Marcq / Mark

Tel : +32(0)2 396 19 10 Fax : +32(0)2 396 06 73

info@aescientific.be

www.aescientific.be

*Price valid until 31-12-2025

All prices are w/o VAT