

EX Pell & EX Pell *Secure*

MICROCENTRIFUGE TUBES

CAPP **EX Pell** and **EX Pell Secure MCT** are manufactured from virgin PP, out of the finest molds available.

All CAPP **EX Pell** and **EX Pell Secure MCT** are **DNase**, **RNase** and **Pyrogen-free** and withstand up to 20.000g.

EX Pell Secure Low-Retention Microcentrifuge Tubes offer a superior smooth inner surface, minimizing sample binding.

EX Pell Secure Black Microcentrifuge Tubes are suitable for applications requiring photo-sensitive, fluorescent or ethidium bromide probes. Unlike amber tinted tubes, they are completely opaque and protect even the most light-sensitive samples.



EX Pell

Microtube
5101500C

EX Pe *Secure*

Microtube
5101500C

Ordering Information 2024

5001500C	ExpellPlus Microcentrifuge tubes 1.5mL, bulk, low retention , 20x500 pcs.	270.00 € / 10 000 pcs
5101500C	ExpellSecure Microcentrifuge tubes 1.5mL, bulk, bag, 20x500 pcs.	172.00 € / 10 000 pcs
5101501C	Expell Microcentrifuge tubes 1.5mL, bulk bag, 20x500 pcs.	172.00 € / 10 000 pcs
5101505C	Expell Secure Microcentrifuge tubes 1.5 mL, pre-sterile, bulk, 20x500 pcs.	182.00 € / 10 000 pcs
5101550C	ExpellSecure Microcentrifuge tubes 1.5mL, bulk, black , 20x500 pcs.	186.00 € / 10 000 pcs
5102000C	ExpellSecure Microcentrifuge tubes 2.0mL, bulk, bag, 16x500 pcs.	182.00 € / 8 000 pcs
5102001C	Expell Microcentrifuge tubes 2.0mL, bulk, bag, 10x500 pcs.	126.00 € / 5 000 pcs
5102005C	ExpellSecure microcentrifuge tubes 2.0 mL, pre-sterile bulk, 16x500 pcs.	194.00 € / 8 000 pcs
5101590C	ExpellSecure Microcentrifuge tubes 1.5mL, bulk, assorted colours, 20x500 pcs. (magenta, blue, yellow, green, orange)	186.00 € / 10 000 pcs

*All prices are without VAT

Made in Germany

Order by fax: 0032 (2) 396 06 73 or email: order@aescientific.be