



CAPP Controller

BUY 2 Capp Controller

Get 1 Case of Serological pipette 5 ml (SP-5-C)

FREE

MOTORIZED PIPETTE CONTROLLER

Features:

- Light and fatigue free: weighs only 167g including the battery
- Lithium battery with no memory effect and up to 12 hours of continuous use
 - Can be charged while in use
- Powerful: 50mL pipette is filled in less than 10 seconds
 - Express dispensing: depress the dispensing button beyond the built-in point of resistance
- Gravity delivery and blow-out dispensing mode selection
 - Fully adjustable speed-setting
- Exchangeable 0.2 μ m PTFE disc filter acts as effective shut-down system if overflow occurs
 - Exchangeable pipette holders are made of silicon
- For use with glass or plastic pipettes ranging from 0.1mL to 100mL
 - Vapour outlet for venting out corrosive liquids or fumes
 - Fully autoclavable nose cone
 - UV-resistant housing



- Complies with CE and RoHS

Ordering Information 2019 - 2020

Cat. No. Description

PA-100 CAPPController Motorized Controller, 0.1-100mL

225.00 €

* Price is without VAT

www.aescientific.com ++32 2 396 19 10 info@aescientific.com