

CAPP *Tempo*

ELECTRONIC PIPETTES

The **CAPP *Tempo*** pipette controller has revolutionary features that separate it from other pipette controllers.

LED Display

- Display shows aspiration and dispense speeds
- Displays remaining battery life and flashes when the power is low
- Display turns-off after 20 seconds of non-use to conserve battery

Adjustable speeds

- Powerful motor fills 25mL in 5 seconds (at speed level six)
- Separate speed controls for aspiration and dispense modes
- Six different speed settings for both aspiration and dispense modes

Electronic Brake System (EBS)

- Apply EBS to slow down at the last second to dial in the most accurate volume
- Press EBS a second time to return to max speed
- EBS maximizes speed and precision

Power Supply

- Lithium ion battery, full charge provides 8 hours of continuous use
- Battery can be charged while pipette is in use

CAPPTempo
&
Charging Stand

Ordering Information 2021

* All prices are without VAT, offer valid only in Belgium & Luxembourg

295.00 €

Cat. No.	Description
T100-B	CAPPTempo Blue
T100-G	CAPPTempo Green
T100-R	CAPPTempo Red
T100-V	CAPPTempo Violet
T100-CS	CAPPTempo Charging Stand

